

Rabbit Anti-Ing3 Polyclonal Antibody For Wb And Elisa - Q9Nxr8

Item Number: CM0024792



Introduction

Rabbit polyclonal antibody against human ING3, validated for WB and ELISA. Reacts with human, mouse, and rat. Detects ING3, a NuA4 complex subunit, supporting chromatin and transcription research.

[Learn More](#)

Parameter	Value
Product Name	ING3 Rabbit pAb
Remarks/Alias	Eaf4; ING2; MEAF4; p47ING3; ING3
Species	Human
GeneID	54556
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human ING3 (NP_938008.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXR8
Protein Weight	47kDa
Shipping	Ice bag