


Rabbit Anti-Kiaa0101 Polyclonal Antibody For Wb And Elisa - Q15004

Item Number: CM0013544



Introduction

Rabbit polyclonal antibody targeting human KIAA0101 (PAF15), validated for Western blot and ELISA. Suitable for DNA repair, cell cycle, and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	KIAA0101 Rabbit pAb
Remarks/Alias	L5, PAF, OEATC, PAF15, OEATC1, p15PAF, NS5ATP9, OEATC-1, p15/PAF, KIAA0101, p15(PAF)
Species	Human
Gene ID (Human)	9768
Gene ID	9768
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 39-111 of human KIAA0101 (NP_055551.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15004
Protein Weight	12kDa
Shipping	Ice bag