


Rabbit Anti-Megalin/Lrp2 Monoclonal Antibody For If-P And Fc (Intra) - P98164

Item Number: CM0009403



Introduction

Rabbit monoclonal antibody targeting human Megalin/LRP2. Validated for immunofluorescence on paraffin sections (IF-P) and intracellular flow cytometry (FC). Cross-reacts with mouse and rat. Suitable for studies of endocytosis, renal reabsorption, and neurodevelopment.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Megalin/LRP2 Monoclonal Antibody
Remarks/Alias	LRP2, DBS, GP330, LDL receptor related protein 2
Species	Human
GeneID (Human)	4036
GeneID	4036
Immunogen	Recombinant fusion protein corresponding to amino acids 4446-4655 of human Megalin/LRP2 (NP_004516.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, FC (intra)
Cross Reactivity	Human, Mouse, Rat
SWISS	P98164
Protein Weight	521kDa
Shipping	Ice bag