

Rabbit Anti-Monkey Ige Monoclonal Antibody (Mab) - Ige

Item Number: CM0005502

Introduction

Rabbit monoclonal antibody against Cynomolgus monkey IgE, produced using recombinant fusion protein immunogen. Suitable for flow cytometry applications in diagnostic and research settings.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Monkey IgE mAb
Remarks/Alias	IgE; IGHE
Species	Cynomolgus monkey
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-429 of Monkey IgE(EHH62331.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cynomolgus
Shipping	Ice bag