

Rabbit Anti-Mouse Cd10 Mab For Fc - Q61391

Item Number: CM0007975



Introduction

-conjugated rabbit monoclonal antibody against mouse CD10 (neprilysin) for flow cytometry. Recombinant immunogen, validated in mouse. Key tool for neuropeptide and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD10 mAb
Remarks/Alias	Mme
Species	Mouse
GeneID (Human)	4311
GeneID	17380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of mouse CD10 (NP_032630.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61391
Protein Weight	86kDa
Shipping	Ice bag