

Rabbit Anti-Mouse Cd10 (Mme) Monoclonal Antibody For Flow Cytometry - Q61391

Item Number: CM0004322



Introduction

-conjugated rabbit monoclonal antibody against mouse CD10 (neprilysin, Mme). Validated for flow cytometry (FC). Ideal for studying neuropeptide degradation and cell surface marker expression in mouse models. - Q61391

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD10 mAb
Remarks/Alias	Mme
Species	Mouse
GeneID (human)	4311
GeneID	17380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of mouse CD10 (NP_032630.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61391
Protein Weight	86kDa
Shipping	Ice bag