

Rabbit Anti-Mouse Cd10 Monoclonal Antibody For Flow Cytometry - Q61391

Item Number: CM0005220

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD10 (neprilysin) for flow cytometry. Validated in mouse samples, recognizes an 86kDa cell surface metalloprotease involved in opioid peptide and ANF degradation. Ideal for immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD10 mAb
Remarks/Alias	Mme
Species	Mouse
GeneID (Human)	4311
GeneID (Mouse)	17380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of mouse CD10 (NP_032630.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61391
Protein Weight	86kDa
Shipping	Ice bag