

Rabbit Anti-Mouse Cd102 Monoclonal Antibody For Flow Cytometry - P35330

Item Number: CM0007926



Introduction

Rabbit monoclonal antibody against mouse CD102 (ICAM-2), conjugated to for flow cytometry. Recognizes mouse CD102, a 31kDa cell adhesion molecule involved in leukocyte and endothelial interactions. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD102 mAb
Remarks/Alias	CD102; Ly-60; Icam-2
Species	Mouse
GeneID (Human)	3384
GeneID (Mouse)	15896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-222 of mouse CD102 (NP_034624.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35330
Protein Weight	31kDa
Shipping	Ice bag