

Rabbit Anti-Mouse Cd102 Monoclonal Antibody For Flow Cytometry - P35330

Item Number: CM0006084



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD102 (Icam-2) for flow cytometry (FC). Reactive with mouse samples; ideal for studying leukocyte adhesion, immune surveillance, and endothelial cell interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD102 mAb
Remarks/Alias	CD102; Ly-60; Icam-2
Species	Mouse
GeneID (human)	3384
GeneID	15896
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-222 of mouse CD102 (NP_034624.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35330
Protein Weight	31kDa
Shipping	Ice bag