

Rabbit Anti-Mouse Cd103 Monoclonal Antibody For Flow Cytometry - Q60677

Item Number: CM0010016

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD103 (integrin alpha-E). Validated for flow cytometry, reacts with mouse. Useful for studying lymphocyte-epithelial adhesion. Molecular weight ~129 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD103 mAb
Remarks/Alias	CD103; aM290; alpha-E1; A530055J10; alpha-M290
Species	Mouse
GenelD (human)	3682
GenelD	16407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 191-385 of mouse CD103 (NP_032425.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60677
Protein Weight	129kDa
Shipping	Ice bag