

Rabbit Anti-Mouse Cd105/Endoglin Mab For Fc - Q63961

Item Number: CM0011954

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD105/Endoglin. Suitable for flow cytometry. Enables study of angiogenesis, TGF-beta/BMP signaling, and vascular biology. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD105/Endoglin mAb
Remarks/Alias	Endo; CD105; S-endoglin
Species	Mouse
Gene ID (Human)	2022
Gene ID	13805
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q63961
Protein Weight	70kDa
Shipping	Ice bag