

Rabbit Anti-Mouse Cd105/Endoglin Mab For Flow Cytometry - Q63961

Item Number: CM0005422



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD105/Endoglin (Q63961). Recombinant immunogen for specificity. Ideal for flow cytometry-based angiogenesis and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD105/Endoglin mAb
Remarks/Alias	Endo; CD105; S-endoglin
Species	Mouse
GenElD (Human)	2022
GenElD	13805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-581 of mouse CD105 (NP_031958.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q63961
Protein Weight	70kDa
Shipping	Ice bag