

Rabbit Anti-Mouse Cd112/Nectin-2 Monoclonal Antibody For Flow Cytometry - P32507

Item Number: CM0008884



Introduction

-conjugated rabbit anti-mouse CD112/Nectin-2 monoclonal antibody for flow cytometry. Targets T-cell modulator and adhesion protein nectin-2; suitable for mouse samples in research on immune regulation and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD112/Nectin-2 mAb
Remarks/Alias	MPH; Pvr; Pvs; Cd112; Pvrl2; nectin-2
Species	Mouse
Gene ID (Human)	5819
Gene ID	19294
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 32-345 of mouse CD112/Nectin-2 (NP_033016.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P32507
Protein Weight	57kDa
Shipping	Ice bag