

Rabbit Anti-Mouse Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P05532

Item Number: CM0011525



Introduction

Flow-cytometry-validated -conjugated rabbit monoclonal antibody targeting mouse CD117/c-Kit (SWISS: P05532). Detects surface receptor on hematopoietic and germ cells; suitable for stem cell research and mast cell analysis.

[Learn More](#)

Specification	Detail
Product Name	Rabbit anti-Mouse CD117/c-Kit mAb
Remarks/Alias	W, Bs, Fdc, Ssm, SCO1, SCO5, SOW3, CD117, c-KIT, Tr-kit, Gsfsc01, Gsfsc05, Gsfscow3
Species	Mouse
GeneID (Human)	3815
GeneID	16590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117/c-Kit (NP_001116205.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	22kDa/108kDa/109kDa
Shipping	Ice bag