

Rabbit Anti-Mouse Cd119/Ifn-Gamma Receptor 1 Mab - P15261

Item Number: CM0011718

Introduction

-conjugated rabbit monoclonal antibody against mouse CD119 (IFN-gamma receptor 1). Suitable for flow cytometry to characterize interferon gamma signaling in immune cells. Specific to mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD119/IFN-gamma receptor 1 mAb
Remarks/Alias	Ifgr; CD119; Ifngr; Nktar; IFN-gammaR
Species	Mouse
GeneID (human)	3459
GeneID	15979
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-250 of mouse CD119/IFN-gamma receptor 1 (NP_034641.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15261
Protein Weight	52 kDa
Shipping	Ice bag