

Rabbit Anti-Mouse Cd11A Monoclonal Antibody For Flow Cytometry - P24063

Item Number: CM0007736

Introduction

This -conjugated rabbit monoclonal antibody specifically recognizes mouse CD11a (LFA-1 alpha chain), validated for flow cytometry (FC). Ideal for immunophenotyping and leukocyte adhesion studies in murine models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD11a mAb
Remarks/Alias	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Species	Mouse
GeneID (Human)	3683
GeneID (Mouse)	16408
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24063
Protein Weight	128kDa
Shipping	Ice bag