


Rabbit Anti-Mouse Cd11B Monoclonal Antibody For Flow Cytometry - P05555

Item Number: CM0008731



Introduction

Rabbit monoclonal antibody targeting mouse CD11b (Integrin alpha-M, CR3, MAC-1). Validated for flow cytometry (FC) applications. Recombinant protein immunogen (aa 143-337). Cross-reacts with mouse. Ideal for monocyte/macrophage and microglia research. UniProt P05555.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik
Species	Mouse
GeneID	16409
GeneID (Human)	3684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	127 kDa
Shipping	Ice bag