

Rabbit Anti-Mouse Cd123 Monoclonal Antibody For Fc - P26952

Item Number: CM0009282

Introduction

-conjugated rabbit monoclonal antibody against mouse CD123 (IL-3 receptor alpha), validated for flow cytometry (FC). Suitable for detecting mouse hematopoietic progenitor cells and IL-3 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD123 mAb
Remarks/Alias	CD123; SUT-1; CDw123
Species	Mouse
Gene ID (human)	3563
Gene ID	16188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-331 of mouse CD123 (NP_032395.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P26952
Protein Weight	43kDa
Shipping	Ice bag