

Rabbit Anti-Mouse Cd124/Il-4Ra Mab For Fc - P16382

Item Number: CM0011653

Introduction

-conjugated rabbit anti-mouse CD124/IL-4Rα monoclonal antibody designed for flow cytometry. Targets the extracellular domain of mouse IL-4R, enabling studies of IL-4/IL-13 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD124/IL-4Ra mAb
Remarks/Alias	Il4r; CD124
Species	Mouse
Gene ID (human)	3566
Gene ID	16190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-233 of mouse CD124/IL-4Ra (NP_001008700.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16382
Protein Weight	87kDa
Shipping	Ice bag