

Rabbit Anti-Mouse Cd124/Il-4Ra Monoclonal Antibody For Fc - P16382

Item Number: CM0010156

Introduction

-conjugated rabbit monoclonal antibody against mouse CD124/IL-4R α for flow cytometry. Recognizes extracellular domain (aa 26-233), validated for mouse samples. Ideal for studying IL-4/IL-13 signaling in allergy and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD124/IL-4Ra mAb
Remarks/Alias	Il4r; CD124
Species	Mouse
GeneID (Human)	3566
GeneID	16190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-233 of mouse CD124/IL-4Ra (NP_001008700.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16382
Protein Weight	87kDa
Shipping	Ice bag