

Rabbit Anti-Mouse Cd126/IL-6R α Chain Mab - P22272

Item Number: CM0009918

Introduction

Recombinant rabbit monoclonal antibody targeting mouse IL-6R α (CD126), conjugated to for flow cytometry. Detects native IL-6R α , supports immunology and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD126/IL-6R α chain mAb
Remarks/Alias	IL6r; CD126; IL-6R; IL-6RA; IL-6R-alpha
Species	Mouse
Gene ID (Human)	3570
Gene ID (Mouse)	16194
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 20-357 of mouse CD126/IL-6R α chain (NP_034689.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P22272
Protein Weight	50kDa
Shipping	Ice bag