

Rabbit Anti-Mouse Cd126/Il-6Rα Chain Monoclonal Antibody For Fc - P22272

Item Number: CM0005512

Introduction

-conjugated rabbit monoclonal antibody against mouse IL-6 receptor alpha (CD126/IL-6Rα). Validated for flow cytometry (FC). Suitable for studying IL-6 signaling, immune response, and inflammatory diseases. UniProt: P22272.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD126/IL-6Rα chain mAb
Remarks/Alias	Il6r; CD126; IL-6R; IL-6RA; IL-6R-alpha
Species	Mouse
GenelD (human)	3570
GenelD (mouse)	16194
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-357 of mouse CD126/IL-6Rα chain (NP_034689.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P22272
Protein Weight	50kDa
Shipping	Ice bag