

Rabbit Anti-Mouse Cd126/Il-6Rα Chain Monoclonal Antibody For Flow Cytometry - P22272

Item Number: CM0009141

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD126/IL-6 receptor alpha (IL-6Rα) for flow cytometry. Detects both membrane-bound and soluble forms, ideal for immune response and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD126/IL-6Rα chain mAb
Remarks/Alias	Il6r; CD126; IL-6R; IL-6RA; IL-6R-alpha
Species	Mouse
GeneID (Human)	3570
GeneID	16194
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-357 of mouse CD126/IL-6Rα chain (NP_034689.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P22272
Protein Weight	50kDa
Shipping	Ice bag