

Rabbit Anti-Mouse Cd127/IL-7R α Monoclonal Antibody For Fc - P16872

Item Number: CM0009366



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD127/IL-7R α , validated for flow cytometry. Specific for mouse, with a predicted molecular weight of 51 kDa. Suitable for immunology and cytokine receptor research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD127/IL-7R α mAb
Remarks/Alias	CD127; IL-7Ralpha
Species	Mouse
Gene ID (Human)	3575
Gene ID (Mouse)	16197
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of mouse CD127/IL-7R α (NP_612564.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16872
Protein Weight	51kDa
Shipping	Ice bag