

Rabbit Anti-Mouse Cd127/IL-7R α Monoclonal Antibody For Flow Cytometry - P16872

Item Number: CM0010831



Introduction

High-quality -conjugated rabbit monoclonal antibody against mouse CD127/IL-7R α , suitable for flow cytometry. Ideal for investigating interleukin-7 and TSLP receptor signaling in spleen, thymus, and immune cell development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD127/IL-7R α mAb
Remarks/Alias	CD127; IL-7R α
Species	Mouse
Gene ID (Human)	3575
Gene ID (Mouse)	16197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of mouse CD127/IL-7R α (NP_612564.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Mouse
SWISS	P16872
Protein Weight	51kDa
Shipping	Ice bag