


Rabbit Anti-Mouse Cd130/Gp130 Monoclonal Antibody For Flow Cytometry - Q00560

Item Number: CM0006259



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD130/gp130 (IL6ST) extracellular domain. Validated for flow cytometry, reactive with mouse samples. Ideal for immune, bone, and neuronal signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD130/gp130 mAb
Remarks/Alias	CD130; gp130; D13Ert699e; 5133400A03Rik
Species	Mouse
Gene ID (Human)	3572
Gene ID (Mouse)	16195
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-617 of mouse CD130/gp130 (NP_034690.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00560
Protein Weight	102kDa
Shipping	Ice bag