

## Rabbit Anti-Mouse Cd138/Syndecan-1 Monoclonal Antibody For If-P And Fc - P18828

Item Number: CM0012444

### Introduction

-conjugated rabbit anti-mouse CD138/Syndecan-1 monoclonal antibody validated for immunofluorescence (IF-P) and flow cytometry (FC). Suitable for studying exosome biogenesis, cell adhesion, and signaling.

[Learn More](#)

| Parameter        | Value                                                                                                                                             |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD138/Syndecan-1 mAb                                                                                                            |
| Remarks/Alias    | Sstn; Synd; CD138; Synd1; syn-1                                                                                                                   |
| Species          | Mouse                                                                                                                                             |
| GeneID (Human)   | 6382                                                                                                                                              |
| GeneID           | 20969                                                                                                                                             |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-255 of mouse CD138/Syndecan-1 (NP_035649.1). |
| Source           | Rabbit                                                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                                                             |
| Application      | IF-P, FC                                                                                                                                          |
| Cross Reactivity | Mouse                                                                                                                                             |
| SWISS            | P18828                                                                                                                                            |
| Protein Weight   | 33kDa                                                                                                                                             |
| Shipping         | Ice bag                                                                                                                                           |