


Rabbit Anti-Mouse Cd14 Monoclonal Antibody For Fc - P10810

Item Number: CM0010592



Introduction

-conjugated rabbit monoclonal antibody against mouse CD14, validated for flow cytometry. Mouse CD14 is a key innate immunity coreceptor for LPS and bacterial lipopeptides. Ideal for immunophenotyping macrophage subsets in immunological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD14 mAb
Remarks/Alias	CD14
Species	Mouse
GeneID (Human)	929
GeneID (Mouse)	12475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-345 of mouse CD14 (NP_033971.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P10810
Protein Weight	39 kDa
Shipping	Ice bag