


Rabbit Anti-Mouse Cd14 Monoclonal Antibody For Fc - P10810

Item Number: CM0010270



Introduction

-conjugated rabbit monoclonal antibody against mouse CD14, validated for flow cytometry (FC). Targets the 39 kDa monocyte differentiation antigen, a co-receptor for bacterial lipopolysaccharide (LPS) and innate immune response mediator.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD14 mAb
Remarks/Alias	CD14
Species	Mouse
GeneID (Human)	929
GeneID	12475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-345 of mouse CD14 (NP_033971.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P10810
Protein Weight	39kDa
Shipping	Ice bag