


Rabbit Anti-Mouse Cd140B/Pdgfr Beta Monoclonal Antibody For Fc - P05622

Item Number: CM0009106



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD140b/PDGFR beta receptor. Validated for flow cytometry applications, reactive specifically with mouse. Predicted molecular weight 123 kDa. Ideal for studies of PDGF signaling in development, angiogenesis, and renal fibrosis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD140b/PDGFR beta mAb
Remarks/Alias	Pdgfr; CD140b; PDGFR-1
Species	Mouse
GeneID (Human)	5159
GeneID (Mouse)	18596
Immunogen	Recombinant protein (recombinant fusion protein containing amino acids 32-530 of mouse CD140b/PDGFR beta, NP_032835.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05622
Protein Weight	123kDa
Shipping	Ice bag