

Rabbit Anti-Mouse Cd144/Ve-Cadherin Monoclonal Antibody For Flow Cytometry - P55284

Item Number: CM0008033



Introduction

ABflo 647-conjugated rabbit monoclonal antibody targeting mouse CD144/VE-Cadherin (SWISS P55284) for flow cytometry. Enables detection of vascular endothelial cadherin in endothelial cell adhesion, junction organization, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD144/VE-Cadherin mAb
Remarks/Alias	7B4; Vec; VECD; Cd144; VEcad; VE-Cad
Species	Mouse
GeneID (Human)	1003
GeneID (Mouse)	12562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 46-592 of mouse CD144/VE-Cadherin (NP_033998.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P55284
Protein Weight	88kDa
Shipping	Ice bag