

Rabbit Anti-Mouse Cd146 Monoclonal Antibody For Fc - Q8R2Y2

Item Number: CM0010547

Introduction

Conjugated rabbit monoclonal antibody against mouse CD146 (MUC18/gicerin). Validated for flow cytometry. Targets the extracellular domain (aa 24-563). Ideal for melanoma adhesion and endothelial monolayer studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD146 mAb
Remarks/Alias	CD146; CD149; Muc18; s-endo; 1-gicerin; s-gicerin
Species	Mouse
GeneID (Human)	4162
GeneID	84004
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-563 of mouse CD146 (NP_075548.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8R2Y2
Protein Weight	71kDa
Shipping	Ice bag