

Rabbit Anti-Mouse Cd161C/Nk1.1 Monoclonal Antibody For Flow Cytometry - P27814

Item Number: CM0011917

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD161c/NK1.1 (Klrb1c), validated for flow cytometry. Recombinant immunogen ensures specific detection of NK cells. Ideal for immunology and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD161c/NK1.1 mAb
Remarks/Alias	Nk1; Ly59; Nk-1; CD161; Ly-59; Ly55c; NK1.1; NKRP1; Nk1.2; Klrb1b; NK-RP1; Nk-1.2; Nkrp1b; Nkrp1c; ly-55c; NKRP140; NKR-P1.9
Species	Mouse
Gene ID (Human)	3820
Gene ID	17059
Immunogen	Recombinant Protein (Recombinant fusion protein containing a sequence corresponding to amino acids 67-223 of mouse CD161c/NK1.1 (NP_001153376.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P27814
Protein Weight	25kDa
Shipping	Ice bag