

Rabbit Anti-Mouse Cd161C/Nk1.1 Monoclonal Antibody For Flow Cytometry - P27814

Item Number: CM0007954

Introduction

-conjugated rabbit monoclonal antibody against mouse CD161c/NK1.1 for flow cytometry. Detects NK cell surface antigen NK1.1; suitable for immunophenotyping and functional studies of natural killer cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD161c/NK1.1 mAb
Remarks/Alias	Nk1, Ly59, Nk-1, CD161, Ly-59, Ly55c, NK1.1, NKRP1, Nk1.2, Klrb1b, NK-RP1, Nk-1.2, Nkrp1b, Nkrp1c, ly-55c, NKRP140, NKR-P1.9
Species	Mouse
Gene ID (Human)	3820
Gene ID (Mouse)	17059
Immunogen	Recombinant fusion protein containing amino acids 67-223 of mouse CD161c/NK1.1 (NP_001153376.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P27814
Protein Weight	25kDa
Shipping	Ice bag