

## Rabbit Anti-Mouse Cd162/Psgl-1 Monoclonal Antibody For Fc - Q3Ta56

Item Number: CM0012778

### Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD162/PSGL-1 (Selp1g). Designed for flow cytometry applications, cross-reacts with mouse. Membrane protein with signal peptide and transmembrane domain.

[Learn More](#)

| Parameter        | Value                              |
|------------------|------------------------------------|
| Product Name     | Rabbit anti-Mouse CD162/PSGL-1 mAb |
| Remarks/Alias    | CD162; Psgl1; Selp1; Selp1; Psgl-1 |
| Species          | Mouse                              |
| Gene ID (Human)  | 6404                               |
| Gene ID          | 20345                              |
| Immunogen        | Recombinant protein                |
| Source           | Rabbit                             |
| Category         | Monoclonal Antibodies              |
| Application      | FC                                 |
| Cross Reactivity | Mouse                              |
| SWISS            | Q3TA56                             |
| Protein Weight   | 44kDa                              |
| Shipping         | Ice bag                            |