

Rabbit Anti-Mouse Cd162/Psgl-1 Monoclonal Antibody For Fc - Q3Ta56

Item Number: CM0011161

Introduction

Fluorophore-conjugated rabbit monoclonal antibody specifically targeting mouse CD162/PSGL-1 (Q3TA56). Validated for flow cytometry (FC), it serves as a reliable tool for immunological studies of cell adhesion and leukocyte trafficking.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD162/PSGL-1 mAb
Remarks/Alias	CD162; Psgl1; Selp1; Selpl; Psgl-1
Species	Mouse
Gene ID (Human)	6404
Gene ID	20345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-320 of mouse CD162/PSGL-1 (NP_033177.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q3TA56
Protein Weight	44kDa
Shipping	Ice bag