

Rabbit Anti-Mouse Cd19 Mab For Fc - P25918

Item Number: CM0004657



Introduction

High-quality -conjugated rabbit monoclonal antibody against mouse CD19, validated for flow cytometry (FC). Enables precise B-cell identification in spleen, bone marrow, and peripheral blood, supporting immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD19 mAb
Species	Mouse
GeneID (Human)	930
GeneID	12478
Immunogen	Recombinant protein of mouse CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag