

Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Fc - P25918

Item Number: CM0010741

Introduction

-conjugated rabbit anti-mouse CD19 monoclonal antibody for flow cytometry (FC). Targets the B-lymphocyte antigen CD19 (60 kDa), a BCR coreceptor critical for B cell development and function. Validated on mouse B cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD19 mAb
Remarks/Alias	Cd19
Species	Mouse
GeneID (Human)	930
GeneID	12478
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-287 of mouse CD19 (NP_033974.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag