

Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Fc - P25918

Item Number: CM0010339

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD19, a B-lymphocyte antigen coreceptor for BCR signaling. Optimized for flow cytometry on mouse samples. Ideal for B-cell immunophenotyping and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD19 mAb
Remarks/Alias	Cd19
Species	Mouse
Gene ID (Human)	930
Gene ID	12478
Immunogen	Recombinant protein of mouse CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag