

Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry - P25918

Item Number: CM0008988

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD19. Validated for flow cytometry, suitable for B-cell immunophenotyping in mouse models. SWISS: P25918.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD19 mAb
Remarks/Alias	CD19
Species	Mouse
GenelD (Human)	930
GenelD	12478
Immunogen	Recombinant protein of mouse CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag