

## Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry - P25918

Item Number: CM0007966

### Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD19 for flow cytometry. Detects B-lymphocyte antigen on splenocytes, bone marrow, and lymph node cells. Recombinant protein-derived, specific for mouse CD19.

[Learn More](#)

| Parameter        | Value                              |
|------------------|------------------------------------|
| Product Name     | Rabbit anti-Mouse CD19 mAb         |
| Species          | Mouse                              |
| GenelD (Human)   | 930                                |
| GenelD           | 12478                              |
| Immunogen        | Recombinant protein of mouse CD19. |
| Source           | Rabbit                             |
| Category         | Monoclonal Antibodies              |
| Application      | FC                                 |
| Cross Reactivity | Mouse                              |
| SWISS            | P25918                             |
| Protein Weight   | 60 kDa                             |
| Shipping         | Ice bag                            |