

Rabbit Anti-Mouse Cd2 Mab For Fc - P08920

Item Number: CM0008183

Introduction

Rabbit anti-Mouse CD2 monoclonal antibody validated for flow cytometry (FC). Specifically detects mouse CD2 (Ly37/LFA-2), suitable for T-cell adhesion and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD2 mAb
Remarks/Alias	Ly37; LFA-2; Ly-37
Species	Mouse
GeneID (Human)	914
GeneID	12481
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-203 of mouse CD2 (NP_038514.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08920
Protein Weight	38kDa
Shipping	Ice bag