

Rabbit Anti-Mouse Cd2 Monoclonal Antibody For Fc - P08920

Item Number: CM0009007

Introduction

-labeled rabbit monoclonal antibody against mouse CD2 (Ly-37/LFA-2) for flow cytometry. Reacts with mouse, immunogen from amino acids 23-203. Ideal for T-cell adhesion and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD2 mAb
Remarks/Alias	Ly37; LFA-2; Ly-37
Species	Mouse
Gene ID (Human)	914
Gene ID	12481
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 23-203 of mouse CD2 (NP_038514.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08920
Protein Weight	38kDa
Shipping	Ice bag