

Rabbit Anti-Mouse Cd2 Monoclonal Antibody For Fc - P08920

Item Number: CM0005714

Introduction

-conjugated rabbit monoclonal antibody against mouse CD2 (LFA-2, Ly37) for flow cytometry. Recognizes the T-cell surface antigen involved in adhesion and activation. Suitable for mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD2 mAb
Remarks/Alias	Ly37; LFA-2; Ly-37
Species	Mouse
GeneID (Human)	914
GeneID	12481
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 23-203 of mouse CD2 (NP_038514.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08920
Protein Weight	38kDa
Shipping	Ice bag