


Rabbit Anti-Mouse Cd200/Ox2 Monoclonal Antibody For Fc - O54901

Item Number: CM0004797



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD200/OX-2 membrane glycoprotein. Validated for flow cytometry, specific to mouse. Unconjugated recombinant immunogen, amino acids 31-232. UniProt O54901.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD200/OX2 mAb
Remarks/Alias	OX2; Mox2
Species	Mouse
GeneID (Human)	4345
GeneID	17470
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-232 of mouse CD200/OX2 (NP_034948.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O54901
Protein Weight	31kDa
Shipping	Ice bag