

Rabbit Anti-Mouse Cd201/Epcr Monoclonal Antibody For Fc - Q64695

Item Number: CM0006239

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD201/EPCR (PROCR). Designed for flow cytometry analysis of the endothelial protein C receptor; immunogen spans amino acids 18-210. Reacts specifically with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD201/EPCR mAb
Remarks/Alias	Ccca; Epcr; Ccd41
Species	Mouse
GeneID (Human)	10544
GeneID	19124
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of mouse CD201/EPCR(NP_035301.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64695
Protein Weight	27kDa
Shipping	Ice bag