

Rabbit Anti-Mouse Cd201/Epcr Monoclonal Antibody For Flow Cytometry - Q64695

Item Number: CM0008628

Introduction

-conjugated rabbit monoclonal antibody against mouse CD201/EPCR (endothelial protein C receptor), validated for flow cytometry. Detects a 27 kDa membrane receptor involved in blood coagulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD201/EPCR mAb
Remarks/Alias	Ccca; Epcr; Ccd41
Species	Mouse
GenelD (Human)	10544
GenelD	19124
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64695
Protein Weight	27kDa
Shipping	Ice bag