

Rabbit Anti-Mouse Cd201/Epcr Monoclonal Antibody For Flow Cytometry - Q64695

Item Number: CM0004004

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD201/EPCR. Validated for flow cytometry (FC), this antibody is ideal for endothelial cell and coagulation research. Molecular weight: 27 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD201/EPCR mAb
Remarks/Alias	Ccca; Epcr; Ccd41
Species	Mouse
Gene ID (Human)	10544
Gene ID	19124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of mouse CD201/EPCR (NP_035301.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64695
Protein Weight	27kDa
Shipping	Ice bag