

Rabbit Anti-Mouse Cd205/Dec-205 Monoclonal Antibody For Fc - Q60767

Item Number: CM0011907

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD205/DEC-205 (Ly75), a 197 kDa endocytic receptor on dendritic cells. Validated for flow cytometry (FC), reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD205/DEC-205 mAb
Remarks/Alias	CD205; DEC205; DEC-205
Species	Mouse
Gene ID (Human)	100526664
Gene ID	17076
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 216-490 of mouse CD205/DEC-205 (NP_038853.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60767
Protein Weight	197kDa
Shipping	Ice bag