

Rabbit Anti-Mouse Cd205/Dec-205 Monoclonal Antibody For Flow Cytometry - Q60767

Item Number: CM0004407

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD205/DEC-205, validated for flow cytometry. Ideal for studying dendritic cell function and antigen uptake pathways in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD205/DEC-205 mAb
Remarks/Alias	CD205; DEC205; DEC-205
Species	Mouse
Gene ID (Human)	100526664
Gene ID	17076
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60767
Protein Weight	197kDa
Shipping	Ice bag